



# AED Maintenance Checklist

<b>AED Location</b>	
<b>Make/ Model</b>	
<b>Serial Number</b>	

	Week	1	2	3	4	5	6	7	8	9	10
<b>AED Condition</b>											
Battery Status/ Date	✓										
Indicator Light Green / Displays "O" for OK	✓										
AED Exterior OK	✓										
<b>AED Components</b>											
Preconnected Pads	✓										
Check expiry date on pads	✓										
CPR Mask	✓										
Clothing Shears	✓										
Safety Razor	✓										
Gloves	✓										
<b>Initial</b>											
<b>Date</b>											



To Order Supplies please visit [eiremed.ie](http://eiremed.ie)